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FABBRICANTE

Euroresearch S.r.l. - Via Larga, 15 - 20122 Milano - Italy
Tel. +39 02 8055660

The medical device Nithya mousse contains heterologous type I collagen. The product takes the form of a foam, packaged in an aluminium container.

Collagen favours the natural proliferation of connective tissue, it induces fibroblast activity and provides the dermis with the scaffolding it requires, thereby recreating the ideal conditions for the collagen fibres to adopt an orderly arrangement.

INDICATIONS

The collagen foam is indicated for use after aesthetic dermatological treatments, such as minor surgery, laser therapy, dermo-abrasion and peeling procedures, once the lesion has closed completely.

INSTRUCTIONS FOR USE

Cleanse the skin before application.

Apply the product 2 to 4 times a day. External use only.

WARNINGS AND PRECAUTIONS

Do not use on broken skin. Do not use on patients with known hypersensitivity to the ingredients of the product.

During pregnancy, the product should be used under medical supervision. If signs of infection occur, contact your doctor, who will prescribe antibiotic therapy if necessary.

Topical tetracycline must not be used in combination with Nithya mousse.

Do not use the product with soaps/detergents and/or disinfectants containing quaternary ammonium salts, silver and other heavy metals, as it could reduce their action. If necessary, ask your doctor or pharmacist which products to use.

Avoid contact with the eyes and mucous membranes. If unintentionally

swallowed or inhaled, contact a doctor immediately. Close the container after each use. Do not use if the packaging is damaged. Do not use the product after the expiry date which is stated on the pack. The expiry date refers to the product in the original packaging, correctly stored. After use, dispose of the container in compliance with applicable national regulations.

KEEP OUT OF THE REACH OF CHILDREN

CONTRAINDICATIONS/ SIDE EFFECTS

Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

STORAGE AND TRANSPORT CONDITIONS

Store the product in a dry place, away from intense heat sources.

Danger: Extremely flammable aerosol. Pressurised container: May burst if heated. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Protect from sunlight. Do not expose to temperatures exceeding 50°C/122°F. Keep out of reach of children. Do not spray on an open flame or other ignition source.

CONTAINER

50ml aluminium can

MANUFACTURER/ DISTRIBUTOR

Euroresearch S.r.l. - Via Larga, 15 - 20122 Milan - Italy
Tel. +39 02 8055660