Nithya Body Cream

COMPOSITION

The medical device Nithya Body Cream contains type I collagen and elastin.

DESCRIPTION

Collagen, together with elastin, helps to provide the necessary support structure to the dermis, moisturizing it and counteracting the relaxation processes, giving tone and elasticity to the skin.

INDICATIONS OF USE

The collagen and elastin cream is indicated for the firming and hydration of all skin types. It is suitable after cosmetic surgery treatments, such as liposuction, abdominoplasty or body sculpture, after complete healing of the skin.

Nithya Body Cream is the perfect solution to our daily agents:

- Barrier fortification: restore & strengthen skin's barrier integrity;
- Anti-dehidratation: improve moisture;
- Dryness reduction: visibly reduce skin dryness aspect;
- Skin relief: soothes skin erythema.

INSTRUCTIONS FOR USE

Apply the product morning and evening, promoting absorption with a gentle massage. For external use only.

WARNINGS AND PRECAUTIONS

Do not use on damaged skin. Do not use on patients with certain individual hypersensitivity to the components of the product. During pregnancy it is recommended to use the product under the supervision

of a doctor. Topical use of tetracycline in combination with Nithya Body Cream should be avoided. Do not use the product with soaps / detergents, disinfectants containing quaternary ammonium salts, silver and other heavy metals, as they could reduce the action.

Avoid contact with eyes and mucous membranes. In case of accidental ingestion or inhalation, seek medical advice immediately. Close the container after each application. Do not use if the package is damaged. Do not use after the expiry date shown on the package. The expiry date refers to an intact and correctly stored package. After use, dispose of the product in accordance to the current national legislation.

KEEP OUT OF REACH OF CHILDREN

CONTRAINDICATIONS / SIDE EFFECTS

No side effects or sensitizations are known or reported.

STORAGE CONDITIONS AND TRANSPORT

To store the product in a dry place away from intense heat.

PACKAGE

Each pack of Nithya Body Cream contains the medical device in a 100 mL tube and instructions for use.

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Euroresearch S.r.l. - Via Larga, 15 - 20122 Milano - Italy